


# PH & ITS MEASUREMENT




**Dr Prasenjit Mitra**

Department of Biochemistry

Sikkim Manipal Institute of Medical Sciences

# Origin of pH

Number of Hydrogen ions ( $H^+$ ) determine acidity or alkalinity




- Number of Hydrogen ions ( $H^+$ ) in water =  $0.0000001 \text{ mol/L}$
- Logarithm of  $H^+$  Concentration
  - $\text{Log } (0.0000001) = \text{log } (10^{-7}) = -7$
- Negative Logarithm of  $H^+$  Concentration
  - $-[\text{Log } (0.0000001)] = -\text{log } (10^{-7}) = -(-7) = 7$



Søren P. L. Sørensen  
(1868-1939)

**p**ower of **H**ydrogen →



pH Scale



# Determination of pH

## Indicators

- Litmus paper
- pH paper

## Colorimeter

## pH meters



# pH meters - History



Arnold Orville Beckman  
(1900-2004)

# pH meters - Types



Handheld pH meter

Bench top pH meter

Continuous in line pH meter





# pH meters – pH Electrode

Glass electrode

Reference electrode

Combination electrode

3 in 1 electrode



# pH meter -- Glass Electrode




Fig Measuring (Glass) electrode

# pH Reference Electrode




Fig Reference electrode

## Reference electrolyte

- Inert
- High ion concentration → Low electrical resistance
- Contact with measuring solution

## Popular Reference Systems

- Mercury/calomel
- Silver/Silver chloride



# pH meter – Combination Electrode




Fig Combination electrode

# pH meter- Working principle


The potential of glass electrode is measured against that of reference electrode

$$E = E^0 + \frac{2.303 RT}{nF} \log a_{H^+}$$

Standard potential when  $a_{H^+} = 1 \text{ mol/L}$

Nernst potential ( $E_N$ )/Slope factor  
Change in potential per pH unit.  
Depends on absolute temperature

Fig Closed circuit of pH meter





# pH Electrode – Working principle




Fig Magnified view of membrane

# pH meter – working principle




Fig Magnified view of membrane with ion model





# Calibration of pH meter

The measuring electrode and reference electrode, when put in a zero solution (7.0 pH buffer) provides a zero mV output.

## Factors causing differences or changes in potential

- Contamination of the reference electrolyte solution.
- Electrolyte evaporation/depletion
- Chemical attack of the silver/silver chloride wire.
- Junction potential.
- Aging of the measuring electrode.

# Calibration of pH meter





# Calibration of pH meter

2 point calibration

Multi point calibration



Fig pH meter with calibrators

# Errors in determination of pH

Alkaline error

Acidic error

Due to reactivity of reference electrolyte

Error due to temperature variation






# Temperature Compensation

Type of Solution	pH value at	
	20°C	30°C
0.001 Mol/L HCl	3.00	3.00
0.001 Mol/L NaOH	11.17	10.83
Phosphate Buffer	7.43	7.40
Tris Buffer	7.84	7.56

Table – Changes in pH with change in temperature

The linear function for temperature versus pH change → 0.003 pH error/pH unit/°C

# Automated temperature compensation (ATC)



# Maintenance & Storage of pH electrode

## Dehydration

Dehydration of  
glass electrode

Dehydration of  
reference  
electrode

## Factors detrimental to electrode life

Chemical attack

Stripping of gel  
layer

## Transport

Avoidance of  
freezing, extreme  
heat, mechanical  
shock and  
vibration

## Storage

At ambient  
temperatures (10-  
30 °C)

Capped

Ideal storage  
solution → 3 -3.5  
M KCl solution



# Definition of pH – a misnomer

- Concentration versus activity
- Activity depends on ionic strength of a solution
- $pH = -\{\log[H^+] \times [f]\}$  where  $f$  is activity co-efficient
- Activity co-efficient depends on total molality of a solution

Molality	0.001	0.005	0.01	0.05	0.1
Activity co-efficient	0.964	0.935	0.915	0.857	0.829

pH of 0.01 M HCl  
=  $-\log(0.01 \times 0.915)$   
= 2.04

pH of 0.01 M HCl with 0.09 M KCl  
=  $-\log(0.01 \times 0.829)$   
= 2.08

pH is negative logarithm of hydrogen ion activity in a solution